
TR

HİPODERMİK İĞNELİ ŞIRINGA KULLANMA TALİMATI

Ürün Tipi

Hipodermik İğneli Şiringa	Hipodermik İğne	
Luer 2 parça Şiringa 2 ml, 2,5 ml, 5 ml, 10 ml, 20ml	14G 15G 16G 18G 19G 20G	
Luer 3 parça Şiringa 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml	21G 22G 23G 24G 25G 26G 27G 29G 30G	
Luer Lock 3 parça Şiringa 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml		

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

Tel: 0212 622 04 00

www.setmedikal.com.tr

HYPODERMIC SYRINGE WITH NEEDLE INSTRUCTION FOR USE

Product Type

Hypodermic Syringe With Needle	Hypodermic Syringe			
Luer 2 parts Syringe 2 ml, 2,5 ml, 5 ml, 10 ml, 20ml	14G	15G	16G	
Luer 3 parts Syringe 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml	18G	19G	20G	
Luer Lock 3 parts Syringe 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml	21G	22G	23G	
	24G	25G	26G	
	27G	29G	30G	

Intended Use

Hypodermic Syringes, made of empty, sterile and disposable plastic materials, are used by end users for the aspiration and injection of liquids after filling. The syringes are primarily for human use. Sterile syringes are designed to be used immediately after filling and are not intended to retain medication for long period of time.

Hypodermic syringes with needle are suitable for use with luer tip connection connectors in accordance with EN ISO 80369-7.

Sterile single use hypodermic syringes with needles are devices used for the introduction and withdrawal of medicinal products and body fluids to and from the body. Syringes are used for intravenous, intramuscular, subcutaneous and intradermal injection. It is also used for the aspiration of the pharmaceutical/drugs from its container in order to administer them to the patient via systems (for ex. via intravenous (IV) port).

Product Components

Name	Raw material
Barrel	PP – Polypropylene
Plunger	PE –Polyethylene, PP – Polypropylene
Gasket	Isoprene rubber, Latex free.
Hypodermic Needle	Stainless Steel, SS304, PP – Polypropylene
Packaging	Medical grade paper and transparent film for EtO (individually blister packed)

Instruction for Use

- Open the package by peeling it from the place marked by the arrow.
- Visually inspect that the content of package is complete and the package is undamaged.
- The product packaging includes the syringe and needle together.
- During priming, hold the syringe with one hand and the needle with the other.
- Make sure that the grading lines on the syringe cylinder are legible.
- Without removing the protective cap of the needle, firmly fix the syringe and needle together by pushing the needle towards the syringe tip and turning it clockwise.
- Carefully remove the needle's protective cap without damaging the needle.
- If liquid is to be injected into the body with a syringe; draw the required amount of liquid into the syringe.
- Make sure the syringe is free of air bubbles. If there is any air bubble in the syringe, be sure to remove it.
- Before injection, check if fluid is flowing in the syringe.
- Inject the liquid to be injected into the vein, intramuscularly, subcutaneously or subcutaneously, in accordance with the application procedure.
- The syringe can also be used for blood aspiration with an appropriately sized hypodermic needle.
- The electronic version of PR.A32 Instruction for use is available at www.setmedikal.com.tr.

Intended Patient Population

There is no limitation for the patient population.

Intended User

Healthcare personnel

Device Description Intended for Use With The Device

Suitable for use with connectors "EN ISO 80369-7 Small diameter connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravenous or hypodermic applications".

Indications

Hypodermic syringes with needle are devices used to administer medical products and liquids. Syringes are designed to be used for intravenous, intramuscular, subcutaneous and intradermal injection, intravenous blood collection and feeding medical systems. The product shelf life is 5 years.

Contraindications

Hypodermic syringes with needle should not be used in more than one patient. The product can be used in a one-time sterile condition with a single patient. Cleaning or resterilization is not possible. It is not intended for use other than its stated intended use. Dispose of as medical waste after use.

Cautions

- Do not use the product if the package is damaged.
- Visually inspect that the content of package is complete and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- Before use, check that the hypodermic needle is not blocked and the flow is not obstructed.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused or re-sterilized.
- Do not touch the needle tip to prevent contamination/ needle sticks. Skin puncture with a contaminated needle may result in serious illness such as hepatitis, AIDS and known/unknown diseases.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- Sterile syringes are designed to be used immediately after filling, so do not store the medicine for a long time.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority

Hypodermic Needle Configurations

Needle Configurations															
Gauge	14G	15G	16G	18G	19G	20G	21G	22G	23G	24G	25G	26G	27G	29G	30G
Color	Pale Green	Blue-Grey	White	Pink	Cream	Yellow	Deep Green	Black	Deep Blue	Purple	Orange	Brown	Medium Grey	Red	Yellow
Needle Length	1"	1 1/4"	1"	1"	1 1/2"	1 1/2"	5/8"	1 1/4"	1"	1"	1"	1/2"	1 1/2"	1/2"	1/2"
	1 1/2"	1 1/2"	1 1/4"	1 1/2"	1"	1"	1 1/2"	1 1/2"	1 1/4"	1 1/4"	5/8"		2"		1/3"

Storage Conditions

- Store at ambient temperature not exceeding 45°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 5 years.

Description of Symbols Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Move vertically
	Catalogue number		Upper limit of temperature (45°C)
	Medical device		Humidity limitation (20% – 80%)
	Unique device identifier		Do not use if package is damaged and consult instruction for use
	Authorized representative in the European Community (United Kingdom)		Do not reuse
	Sterilized using ethylene oxide		Latex free
	Do not resterilize		Non-pyrogenic
	Single sterile barrier system		Professional Use Only
	Single sterile barrier system with protective packaging outside		Consult instruction for use or consult electronic instruction for use
	Caution		
	Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number		

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