


TR

HİPODERMİK İĞNESİZ ŞIRINGA KULLANMA TALİMATI

Ürün Tipi

Hipodermik İğnesiz Şiringa	
Luer 2 parça Şiringa 2 ml, 2,5 ml, 5 ml, 10 ml, 20 ml	
Luer 3 parça Şiringa 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml	
Luer Lock 3 parça Şiringa 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml	

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
Tel: 0212 622 04 00

www.setmedikal.com.tr

HYPODERMICSYRINGES WITHOUT NEEDLES INSTRUCTION FOR USE

Product Type

Hypodermic Syringes Without Needles	
Luer 2 part Syringe 2 ml, 2,5 ml, 5 ml, 10 ml, 20 ml	
Luer 3 part Syringe 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml	
Luer Lock 3 part Syringe 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml	

Intended Use

Sterile single use hypodermic syringes without needles, made of empty, sterile and single use plastic materials, are used by end users for the aspiration and injection of liquids after filling. The syringes are primarily for human use. Sterile syringes are designed to be used immediately after filling and are not intended to retain medication for long period of time.

Hypodermic syringes without needle are suitable for use with luer tip connection connectors in accordance with EN ISO 80369-7.

Sterile single use hypodermic syringes without needles are devices used for the introduction and withdrawal of medicinal products and body fluids to and from the body. If syringes are used in combination with hypodermic needles, they are used for intravenous, intramuscular, subcutaneous and intradermal injection. It is also used for the aspiration of the pharmaceutical/drugs from its container in order to administer them to the patient via systems (for ex. via intravenous (IV) port).

Product Components

Component Name	Raw Material
Barrel	PP – Polypropylene
Plunger	PE – Polyethylene, PP – Polypropylene
Plunger Stopper (Gasket)	Isoprene rubber, Latex free
Packing Material	Medical grade paper and transparent film for EO (individually blister packed)

Instruction for Use

- Open the packaging by stripping it from the opening place indicated by the arrow sign.
- Visually inspect that the content of package is complete and the package is undamaged.
- The product packaging consists of a sterile syringe.
- Make sure that the grading lines on the syringe barrel are legible.
- If liquid is to be injected into the body with a hypodermic needle of the appropriate size with the syringe; withdraw the liquid into the syringe in the required amount.
- Make sure that there are no air bubbles in the syringe; if any, remove them.
- Before injection, check whether liquid flow can be ensured in the syringe.
- Inject the liquid to be injected into the intravenously, intramuscularly, subcutaneously or intradermally, according to the application procedure.
- In case of need, the syringe can also be used for aspiration of blood from the vein with a hypodermic needle of appropriate size.
- The electronic version of PR.A99 Instruction for use is available at www.setmedikal.com.tr.

Intended Patient Population

There is no limitation for the patient population.

Intended User

Healthcare personnel

Description of the device intended for use with the device

It is suitable for use with connection connectors according to "EN ISO 80369-7 Small – bore connectors for liquids and gases in healthcare applications – Part 7: Connectors for intravascular or hypodermic applications"

Indication

Sterile single use hypodermic syringes without needles are devices used for the introduction and withdrawal of medicinal products and body fluids to and from the body. Syringes are used for intravenous, intramuscular, subcutaneous and intradermal injection, and also to feed medical systems. The shelf life of product is 5 years.

Contraindication

Syringes without Hypodermic needles are not used in more than one patient. The product can be used in a one-time sterile condition with a single patient. Cleaning or re-sterilization cannot be performed. It is not intended for use outside the specified intended use. After usage, dispose of the product as medical waste.

Cautions

- Do not use the Product if the package is opened and damaged.
- Visually inspect that the content of package is complete and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused, rework or re-sterilized.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.
- Sterile syringes are designed to be used immediately after filling, so do not store the medicine for a long time.

Storage Conditions

- Store at ambient temperature not exceeding 45°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 5 years.

Description Of The Symbol Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Move vertically
	Catalogue number		Upper limit of temperature (45°C)
	Medical device		Humidity limitation (20% – 80%)
	Unique device identifier		Do not use if package is damaged and consult instruction for use
	Authorized representative in the European Community (United Kingdom)		Do not reuse
	Sterilized using ethylene oxide		Latex free
	Do not resterilize		Non-pyrogenic
	Single sterile barrier system		Professional Use Only
	Single sterile barrier system with protective packaging outside		Consult instruction for use or consult electronic instruction for use
	Caution		
	Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number		

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