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Tel: 0212 622 04 00

www.setmedikal.com.tr

EN Single Use Oral/Enteral Syringe**Product Type****Enteral Syringe Single Use:**





1 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml

Enteral Syringe Single Use Low Dose:

1 ml, 3 ml, 5 ml

Oral Syringe Single Use:

1 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml

Product Type	
Enteral Syringe Single Use Central:	
Enteral Syringe Single Use Low Dose:	
Enteral Syringe Single Use Eccentric:	
Oral Syringe Single Use:	

Intended Use

ENFit Connector Set Medical Oral/Enteral Syringes and Oral Tip Set Medical Oral/Enteral Syringes are designed for the distribution, measurement, and transfer of liquids. They facilitate the delivery of fluids to a patient's gastrointestinal system either orally or enterally. These syringes are intended for use in both clinical and home care settings by users of all age groups, including both healthcare professionals and inexperienced users.

Product Components

Name	Raw material
Barrel	PP – Polypropylene
Plunger	PE – Polyethylene, PP – Polypropylene
Gasket	Isoprene rubber, Latex free
Packaging Material	Medical paper and transparent film suitable for EO sterilization (single unit package)

Instruction for Use

1. For sterile syringes, use aseptic technique. Ensure the plunger is fully depressed prior to filling syringe. For optimal dosage accuracy, attach an feeding connector designed to draw fluid. Fill to prescribed dosage.
2. For accurate measurement, fluid should be filled to black leading edge of the syringe plunger. If using an Feeding Connector, remove and discard. If using a 1 ml, 3 ml or 5 ml Low Dose ENFit® tip syringe, tap out any extra fluid congregating in the tip for improved dosing accuracy. Measurement accuracy is $\pm 5\%$ for volumes up to and including 3 ml, and $\pm 4\%$ for volumes above 3 ml.
3. Attach a tip cap (available separately) after filling, as necessary.
4. Label syringe with patient's name and syringe contents, or follow facility policy on labeling oral/enteral syringes.
5. Remove tip cap, where applicable, and discard before use.
6. Administer fluids orally or by connecting to the appropriate enteral access port, following facility's protocol. Confirm the adapter is connected to an enteral port only and NOT to an IV set.
7. Ensure tight connection to enteral access port prior to use.
8. Dispose syringe in appropriate container after use, following facility protocols.

Intended Patient Population

There are no limitations set for the patient population.

Intended User

Laypersons and healthcare professionals.

Device Description Intended for Use With The Device

Enteral syringes: Suitable for use with connectors "EN ISO 80369-3 Small diameter connectors for liquids and gases in healthcare applications - Part 3: Connectors for enteral applications".

Indications

It is used in cases where the administration of nutrition or medication via the enteral or oral route is necessary. The main indications suitable for its use include:

- Need for enteral nutrition in patients who are unable to take food orally or have inadequate oral intake,
- Individuals with swallowing disorders (dysphagia),
- Patients with impaired consciousness,
- Accurate dosing in cases where oral medication administration is required,
- Individuals unable to feed due to neurological conditions,
- Feeding applications via nasogastric, orogastric, or gastrostomy tubes,
- Enteral medication administration in patients unable to take drugs orally,
- Need for enteral support in neonates and pediatric patients,
- Nutritional support in premature and low birth weight infants,
- Administration of nutritional supplements or support products,
- Initiation of enteral nutrition in the early postoperative period in patients with adequate gastrointestinal function.

Contraindications

The use of enteral and oral syringes is not appropriate in the following situations:

- Conditions requiring parenteral (intravenous) administration,
- Clinical scenarios where the gastrointestinal system is non-functional or severely impaired,

- Cases where oral administration of the drug is contraindicated,
- Patients at high risk of aspiration,
- Use of medications that are not suitable for enteral or oral administration.

⚠ Cautions

- If using a 1 ml, 3 ml or 5 ml Low Dose ENFit® tip syringe, tap out any extra fluid around the moat of the syringe tip for improved dosing accuracy.
- For Oral/Enteral feeding use. ENFit Enteral Syringe with Connector and Oral Tip is designed for enteral fluids only but has the potential to be inadvertently connected to small bore connectors from other healthcare applications. This syringe should not be used with connectors from other healthcare applications.
- If used on a syringe pump, consult pump manufacturer for compatibility and proper instructions for use.
- Re-sterilization and/or re-use may lead to infection or other illness or injury.
- Sterile syringes are sterilized with ethylene oxide gas.
- Do not use syringe if syringe tip is visibly blocked
- Do not use if package is opened or damaged.
- Free from Latex, DEHP, BPA
- After usage, dispose of the Product as medical waste
- The shelf life of product is 5 years.
- Do not use the product after the expiration date printed on the primary packaging.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.
- This syringe is for enteral use only and is suitable for ordinary users, including both healthcare professionals and caregivers. It is designed for tasks such as enteral feeding, medication administration, or flushing feeding tubes.
- Before using the syringe, make sure you understand the correct technique for safe and effective use.
- Check for any fluid leakage from the plunger or connection point, and ensure that the feeding tube is not clogged.
- If you are unsure about how to use the syringe, whether it is compatible with your feeding tube, extension set, or formula, or how and when to clean or dispose of it safely, consult a healthcare professional.
- Incorrect use can lead to feeding errors, aspiration, tube blockages, or contamination. Always follow local guidelines for disposal.
- The product may be used in sterile condition, one-time with one single patient.
- Cleaning or re-sterilization are excluded. Re-using single use products bears potential risks for patient and user.
- Impairment of material properties and inadequate sterility associated with such re-uses may lead to illness, damage or death of the patient.
- Do not use this device in environments that do not meet standard hygiene or temperature conditions (e.g., outdoors, vehicles, or non-clinical settings) unless specifically instructed by a healthcare provider.
- Use only in clean, dry environments below 45°C. Do not expose to direct sunlight or excessive moisture.

Storage Conditions

- Store at ambient temperature not exceeding 45°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 5 years.

Description of Symbols Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Move vertically
	Catalogue number		Upper limit of temperature (45°C)
	Medical device		Humidity limitation (20% – 80%)
	Unique device identifier		Do not use if package is damaged and consult instruction for use

	Sterilized using ethylene oxide		Do not reuse
	Do not resterilize		Latex free
	Single sterile barrier system		Non-pyrogenic
	Single sterile barrier system with protective packaging outside		Consult instruction for use or consult electronic instruction for use
	Caution		Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number
	Marking to meet UK Medical Device Regulation requirements. The number displayed below is the Notified Body number.		

Üretim Yeri

Set Medikal Sanayi ve Ticaret A.Ş.
Osmangazi Mah. Mareşal Fevzi Çakmak Cad. No:18 Esenyurt – İstanbul, Türkiye
Tel: 0212 622 04 00
www.setmedikal.com.tr