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

Product Type	
ENFit Connector – Centric 3 part	
ENFit Connector – Eccentric 3 part	

Intended Use

The reusable enteral feeding syringe with ENFit connector is used to transfer nutrition, medication or water orally, or enterally into the gastrointestinal system of a patient via an enteral feeding tube. The reusable feeding syringes are intended to be used in clinical or home care settings by users ranging from laypersons (under the supervision of a clinician) to clinicians. It is intended to be used multiple times in non-clinical settings by users ranging from clinicians to laypersons in all age groups. The device is indicated for single patient use only.

Product Components

Name	Raw material
Barrel	PP – Polypropylene
Plunger	PE – Polyethylene, PP – Polypropylene
Gasket	Silicone
Packaging	Medical grade paper and transparent film

Instruction for Use

- Ensure the plunger is fully depressed prior to filling syringe. For optimal dosage accuracy, attach an ENFit Connector designed to draw fluid. Fill to prescribed dosage.
- For accurate measurement, fluid should be filled to black leading edge of the syringe plunger. If using an ENFit Connector, remove and discard. If using a 1 ml, 3 ml or 5 ml Low Dose tip syringe, tap out any extra fluid congregating in the tip for improved dosing accuracy. Measurement accuracy is $\pm 5\%$ for volumes up to and including 3 ml, and $\pm 4\%$ for volumes above 3 ml.
- Attach a tip cap (available separately) after filling, as necessary.
- Label syringe with patient's name and syringe contents, or follow facility policy on labeling oral/enteral syringes.
- Remove tip cap, where applicable, and discard before use.
- Administer fluids orally or by connecting to the appropriate enteral access port, following facility's protocol. Confirm the adapter is connected to an enteral port only and NOT to an IV set.
- Ensure tight connection to enteral access port prior to use.

Cleaning Instructions:

Thoroughly wash syringe immediately after each use.

- While syringe tip is submerged in warm soapy (dish soap) water, pull/push plunger the full length of the barrel at least 3 times.
- Remove plunger from barrel. Using your hands, wash the plunger and barrel separately with warm (43-49°C; 110-120°F) soapy water until visible soil is removed. If visible soil remains, repeat cleaning steps using a soft-bristled brush to aid in washing the plunger and the INSIDE of the syringe barrel. Do not use abrasive products on the OUTSIDE of the barrel as this may lead to premature removal of printed markings on the syringe.
- Thoroughly rinse (until no bubbles remain) and allow both parts to dry in a clean area before putting plunger back in barrel. Repeat cleaning steps if visible soil remains.
- Store in a clean, dry area.

Intended Patient Population

There is no limitation for the patient population.

Intended User

Healthcare professionals and laypersons.

Device Description Intended for Use With The Device

Suitable for use with connectors "EN ISO 80369-3 Small diameter connectors for liquids and gases in healthcare applications - Part 3: Connectors for enteral applications".

Indications

Reusable enteral feeding syringes with ENFit connectors are indicated for the repeated administration of nutrition, fluids, medications, or water into the gastrointestinal tract, either orally or via an enteral feeding tube, in a single patient. The device is intended for use in clinical environments (e.g., hospitals, nursing homes) as well as in non-clinical environments (e.g., home care), by healthcare professionals or by caregivers under medical supervision

Contraindications

The device is not suitable for parenteral applications (e.g., intravenous, intramuscular). The syringe must not be used for the administration of blood or bodily fluids. Shared use between multiple patients is strictly contraindicated. Use of syringes presenting mechanical damage (e.g., cracks, deformation, leakage) is contraindicated. The device must not be used in cases where enteral feeding is deemed contraindicated by the physician (e.g., bowel obstruction, risk of perforation, severe gastrointestinal bleeding)

Cautions

- Not following Cleaning Instructions may result in patient harm.
- If using a 1 ml, 3 ml or 5 ml Low Dose tip syringe, tap out any extra fluid around the moat of the syringe tip for improved dosing accuracy.
- For oral/enteral use only. For use with ENFit Connection System. The Enteral Syringe with ENFit Connector is intended for enteral fluids only, but has the potential to misconnect with small bore connectors of other healthcare applications. This syringe should not be used with connectors from other healthcare applications.
- If used on a syringe pump, consult pump manufacturer for compatibility and proper instructions for use.
- Do not use abrasive products on the outside of the syringe barrel as this may lead to premature removal of printed markings on the syringe.
- Do not use syringe if syringe tip is visibly blocked
- Do not use if package is opened or damaged.
- Free from Latex, DEHP, BPA
- After usage, dispose of the product as medical waste
- The shelf life of product is 5 years.
- Do not use the product after the expiration date printed on the primary packaging.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.

- This syringe can be used by both healthcare professionals and ordinary users. If you are unsure about how to use the syringe, how to clean it, whether it is compatible with your diet or formula, or when and how to safely dispose of the syringe, consult a healthcare professional.
- This syringe is for enteral use only and is suitable for ordinary users, including both healthcare professionals and caregivers. It is designed for tasks such as enteral feeding, medication administration, or flushing feeding tubes.
- Before using the syringe, make sure you understand the correct technique for safe and effective use.
- Check for any fluid leakage from the plunger or connection point, and ensure that the feeding tube is not clogged.
- If you are unsure about how to use the syringe, whether it is compatible with your feeding tube, extension set, or formula, or how and when to clean or dispose of it safely, consult a healthcare professional.
- Incorrect use can lead to feeding errors, aspiration, tube blockages, or contamination. Always follow local guidelines for disposal.
- Do not use this device in environments that do not meet standard hygiene or temperature conditions (e.g., outdoors, vehicles, or non-clinical settings) unless specifically instructed by a healthcare provider.
- Use only in clean, dry environments below 45°C. Do not expose to direct sunlight or excessive moisture.
- **Discard syringe in regular household waste when:**
 - Black o-ring gasket is distorted or damaged
 - Barrel markings are difficult to read
 - Plunger is too difficult to move
 - Syringe use exceeds 20 use cycles

Storage Conditions

- Store at ambient temperature not exceeding 45°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 5 years.

Description of Symbols Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Move vertically
	Catalogue number		Upper limit of temperature (45°C)
	Medical device		Humidity limitation (20% – 80%)
	Unique device identifier		Do not use if package is damaged and consult instruction for use
	Non- sterile		Latex free
	Single patient multiple use		Non-pyrogenic
	Caution		Consult instruction for use or consult electronic instruction for use

 2797	Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number	 0086	Marking for meeting the requirements of UK Medical Device Regulation The number displayed on the bottom is the Notified Body Number
---	--	---	--

Manufacturing Site

Set Medikal Sanayi ve Ticaret A.Ş.

Osmangazi Mah. Mareşal Fevzi Çakmak Cad. No:18 Esenyurt – İstanbul Türkiye

Tel: 0212 622 04 00

www.setmedikal.com.tr