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
GB SYRINGE FOR POWER DRIVEN PUMP INSTRUCTION FOR USE**Product Type**

Syringe For Perfusion Pump	Hypodermic Needle
1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml	14G 15G 16G 18G

Product Type	Syringe For Power Driven Pump – Regular Needle
Luer Lock	

Intended Use

Syringes for perfusion pump are intended to be used for injection of medicinal liquids to human body with power driven pumps. If needle is provided, it is used to aspirate medicinal liquid into the syringe.

Syringes for perfusion pump are suitable for use with luer tip connection connectors in accordance with EN ISO 80369-7 and perfusion pumps .

Product Components

Name	Raw material
Barrel	PP – Polypropylene
Plunger	PP – Polypropylene
Gasket	Isoprene rubber, Latex free.
Hypodermic Needle	Stainless Steel (SS304), PP – Polypropylene
Packaging	Medical grade paper and transparent film for EO (individually blister packed)

**Instruction for Use**

- Open the package by peeling it from the place marked by the arrow.
- Visually inspect that the content of package is complete and the package is undamaged.
- During preparation, hold the syringe with one hand and the needle, if available, with the other hand.
- Make sure that the graduation lines on the syringe cylinder are legible.
- Without removing the protective cap of the needle, firmly fix the syringe and needle together by pushing the needle towards the syringe tip and turning it clockwise.
- Carefully remove the needle's protective cap without damaging the needle.
- Draw the required amount of liquid into the syringe.
- Attach the extension line to the syringe and make sure the syringe is free of air bubbles. If there is any air bubble in the syringe, make sure to remove it.
- Before connecting the syringe with perfusion pump, check if fluid is flowing in the syringe and the extension line. Remove air bubbles inside the extension line.
- Connect the syringe to the perfusion pump and use the pump according to the manufacturer's instructions.
- The electronic version of PR.A100 Instruction for use is available at www.setmedikal.com.tr.
- Determine the amount of liquid to be used using the graduated scale on the syringe.

- The capacity tolerance of the graduated scale on the syringe is stated in the table below in accordance with EN ISO 7886-1 Standard.

Nominal capacity of syringe V*, (ml)	Graduated Capacity Tolerance (%)
5≤V≤60	±4

V*: Volume

Intended Patient Population

There is no limitation for the patient population.

Intended User

Healthcare personnel

Device Description Intended for Use With The Device

Suitable for use with connectors "EN ISO 80369-7 Small diameter connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravenous or hypodermic applications".

Perfusion Pump

Indications

Syringes for perfusion pumps are used for injection of medicinal fluids. These syringes are used with power driven pumps. The product shelf life is 5 years.

Contraindications

Syringes for perfusion pump cannot be used in more than one patient. The product can be used in a one-time sterile condition with a single patient. Cleaning or resterilization is not possible. It is not intended for use other than its stated intended use. Dispose of as medical waste after use.

Cautions

- Do not use the product if the package is damaged.
- Visually inspect that the content of package is complete and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- When the sterile barrier system (unit packaging) is opened, connection is made without waiting in accordance with aseptic rules.
- The maximum usage time of the syringe is 24 hours.
- If there is a special warning for the drug administered by syringe, it should be taken into account during application.
- In applications, the capacity tolerance table of the graduated scale given in the instructions for use should be taken into account.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused, rework or re-sterilized.
- Do not touch the needle tip to prevent contamination/needle sticks. Skin puncture with a contaminated needle may result in serious illness such as hepatitis, AIDS and known/unknown diseases.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority

Hypodermic Needle Configurations

Gauge	14G	15G	16G	18G
Color	Pale Green	Blue-Grey	White	Pink
Needle Length	1"	1 1/4"	1"	1"
	1 1/2"	1 1/2"	1 1/4" 1 1/2"	1 1/2"

Storage Conditions

- Store at ambient temperature not exceeding 45°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 5 years.

Description of Symbols Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Move vertically
	Catalogue number		Upper limit of temperature (45°C)
	Medical device		Humidity limitation (20% – 80%)
	Unique device identifier		Do not use if package is damaged and consult instruction for use
	Responsible Person in the United Kingdom		Do not reuse
	Sterilized using ethylene oxide		Latex free
	Do not resterilize		Non-pyrogenic
	Single sterile barrier system		Professional Use Only
	Single sterile barrier system with protective packaging outside		Consult instruction for use or consult electronic instruction for use
	Caution		Marking for meeting the requirements of UK Medical Device Regulation The number displayed on the bottom is the Notified Body Number
	Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number		

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