TR HİPODERMİK İĞNE KULLANIM TALİMATI

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
HYPODERMIC NEEDLE INSTRUCTION FOR USE
Product Type

Hypodermic Needle					
14G	15G	16G	18G	19G	
20G	21G	22G	23G	25G	
26G	27G	29G	30G		

Intended Use

Hypodermic needles made of empty, sterile and disposable stainless steel are used by end users for the injection or withdrawal of both body fluids (blood, etc.) and medical fluids such as medicines. It may also be used for liquid/drug preparation. The needles are primarily for use in humans. There are two types: One is designed with cap and second one has safety cap.

Hypodermic needles are suitable for use with luer tip connection connectors in accordance with EN ISO 80369-7.


With the safety hypodermic needle product, it is aimed to protect medical personnel against possible needle injuries and infection risks.

Product Components

Component Name	Raw Material
Hypodermic Needle	Stainless Steel, SS304, PP - Polypropylene
Protective/Safety Cap	PP – Polypropylene
Packing Material	Medical grade paper and transparent film for EO (individually blister packed)

Instruction for Use

- Open the packaging by stripping it from the opening place indicated by the arrow sign.
- Visually inspect that the content of package is complete and the package is undamaged.
- The product packaging consists of a sterile needle.
- Without removing the protective cap of the needle, firmly secure the syringe and needle together, pushing the needle toward the syringe tip and turning it clockwise.
- Carefully remove the protective cap of the needle without damaging the needle. For safety hypodermic needles, carefully remove the safety cap of the needle without damaging the needle.
- Inject the liquid to be injected into the intravenously, intramuscularly, subcutaneously or intradermally, according to the application procedure.

Hypodermic Safety Needle Usage Image

For hypodermic safety needle, the safety mechanism of hypodermic needles is activated by easily inserting the needle tip removed from the patient's body after the injection into the safety needle cap. When the user hears a click, s/he understands that the safety guard is securely fixed around the needle and is used without any risk of needle sticking. If it is necessary to move the filled syringe to the point of administration, a safe and passive reclosing technique is used to close the needle. Close the safety cap using one of the following methods immediately after the standard procedure:

Push forward with one hand to close the safety cap.

Press your finger on the back of the safety cap (Figure 1.1. – Figure 1.2.) or on a flat surface with the cap to look downwards (Figure 2.) to close the needle and lock the cap.

The safety hypodermic needle is locked when the needle tip is completely closed.

To insert or remove the needle, hold only the colored needle cover (Figure 3.)

- In case of need, hypodermic needles can also be used for aspiration and fluid/drug preparation.
- The electronic version of PR.A78 Instruction for use is available at www.setmedikal.com.tr.

Intended Patient Population

There is no limitation for the patient population.

Intended User

Healthcare personnel

Description of the device intended for use with the device

It is suitable for use with connection connectors according to "EN ISO 80369-7 Small – bore connectors for liquids and gases in healthcare applications – Part 7: Connectors for intravascular or hypodermic applications"

Indication

Hypodermic needles made of empty, sterile and single use stainless steel are used by end users for the injection or withdrawal of both body fluids (blood, etc.) and medical fluids such as medicines. It may also be used for liquid/drug preparation. The shelf life of product is 5 years.

Contraindication

Hypodermic needles are not used in more than one patient. The product can be used in a one-time sterile condition with a single patient. Cleaning or re-sterilization cannot be performed. It is not intended for use outside the specified intended use. After usage, dispose of the product as medical waste.

Cautions

- Do not use the product if the package is opened and damaged.
- Visually inspect that the content of package is complete and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- Before use, check that the hypodermic needle is not blocked and the flow is not obstructed.
- For safety hypodermic needles, activate it away from yourself and others with the one-hand technique. For a high level of safety, use only the finger pad area on the large rough tissue when activating the safety cap. If liquid remains on the needle during the activation of the safety cap or after the injection, it may cause a small amount of splash.
- For safety hypodermic needles, to help avoid HIV (AIDS), HBV (Hepatitis) and other infectious diseases due to accidental needlesticks, activate the protective mechanism immediately after use.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused, rework or re-sterilized.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority

Hypodermic Needle Configurations

Needle Configurations															
Gauge	14G	15G	16G	18G	19G	20G	21G	22G	23G	24G	25G	26G	27G	29G	30G
Color	Pale Green	Blue-Grey	White	Pink	Cream	Yellow	Deep Green	Black	Deep Blue	Purple	Orange	Brown	Medium Grey	Red	Yellow
Needle Length	1"	1 1/4"	1"	1"	1 1/2"	1 1/2"	5/8"	1 1/4"	1"	1"	1"	1/2"	1 1/2"	1/2"	1/2"
	1 1/2"	1 1/2"	1 1/4" 1 1/2"	1 1 1/2"	1"	1"	1 1/2" 1"	1 1/2" 1"	1 1/4"	1 1/4"	5/8"		2" 1/2"		1/3"

Storage Conditions

- Store at ambient temperature not exceeding 45°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 5 years.

Description Of The Symbol Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Move vertically
	Catalogue number		Upper limit of temperature (45°C)
	Medical device		Humidity limitation (20% – 80%)
	Unique device identifier		Do not use if package is damaged and consult instruction for use
	Responsible Person in the United Kingdom		Do not reuse
	Sterilized using ethylene oxide		Latex free
	Do not resterilize		Non-pyrogenic
	Single sterile barrier system		Professional Use Only
	Single sterile barrier system with protective packaging outside		Consult instruction for use or consult electronic instruction for use
	Caution		Marking for meeting the requirements of UK Medical Device Regulation The number displayed on the bottom is the Approved Body Number
	Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number		

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