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Tel: 0212 622 04 00

www.setmedikal.com.tr

BULK HYPODERMIC SYRINGE INSTRUCTION FOR USE
Product Type

Bulk Hypodermic Syringe
Luer 2 part Syringe 2 ml, 2,5 ml, 5 ml, 10 ml, 20 ml
Luer 3 part Syringe 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml
Luer Lock 3 part Syringe 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml

Syringe Type	Bulk Hypodermic Syringe
Luer – Eccentric 2 parts	
Luer – Eccentric 3 parts	
Luer – Centric 2 parts	
Luer – Centric 3 parts	
Luer Lock	

Intended Use

Bulk hypodermic syringes, made of empty, non sterile and single use plastic materials, are generally used to administer or withdraw fluids, with the intended use determined by the end user. Sterilization is performed by the final manufacturer. The product is suitable for EO sterilization. The following parameters must be used for the sterilization process. The final manufacturer is responsible for performing the sterilization process with these validated parameters.

Sterilization Parameters		
EO	Temperature	Humidity
%90 EO, %10 CO ₂	54 – 58°C	40 – 80 %rh

Bulk hypodermic syringes are suitable for use with luer tip connection connectors in accordance with EN ISO 80369-7.

Single use bulk hypodermic syringes are devices used for the introduction and withdrawal of medicinal products and body fluids to and from the body.

Product Components

Component Name	Raw Material
Barrel	PP – Polypropylene
Plunger	PE – Polyethylene, PP – Polypropylene
Plunger Stopper (Gasket)	Isoprene rubber, Latex free
Packing Material	Double PE bag – transport carton

Instruction for Use

The following instructions for use are intended for general use of the syringes, independent of the intended use specified by the final manufacturer.

- Visually check that the contents of the transport carton are complete and that the carton/ bag is undamaged.
- Make sure that the grading lines on the syringe barrel are legible.
- Withdraw the liquid into the syringe in the required amount.
- Before use, check whether liquid flow can be ensured in the syringe.
- Determine the amount of liquid to be used using the graduated scale on the syringe.
- The capacity tolerance of the graduated scale on the syringe is stated in the table below in accordance with EN ISO 7886-1 Standard.

Nominal capacity of syringe V*, (ml)	Graduated Capacity Tolerance (%)
$1 \leq V \leq 3$	± 5
$5 \leq V \leq 60$	± 4

V*: Volume

Intended Patient Population

The patient population for the intended use of the finished product is determined by the final manufacturer.

Intended User

Healthcare personnel

Description of the device intended for use with the device

It is suitable for use with connection connectors according to “EN ISO 80369-7 Small – bore connectors for liquids and gases in healthcare applications – Part 7: Connectors for intravascular or hypodermic applications”

Indication

Non sterile single use bulk hypodermic syringes are devices used for the introduction and withdrawal of medicinal products and body fluids to and from the body. The shelf life of product is 5 years.

Contraindication

It is not intended for use other than that specified by the final manufacturer.
After usage, dispose of the product as medical waste.

Cautions

- Check the product if the transport carton is damaged.
- Do not use the product after the expiration date printed on label.
- If there is a special warning for the drug administered by syringe, it should be taken into account during application.
- In applications, the capacity tolerance table of the graduated scale given in the instructions for use should be taken into account.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused, rework or re-sterilized.
- Cleaning and disinfection cannot be applied to the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.
- The products are suitable for sterilization with EO sterilization method.

Storage Conditions

- Store at ambient temperature not exceeding 45°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 5 years.

Description Of The Symbol Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Move vertically
	Catalogue number		Upper limit of temperature (45°C)
	Medical device		Humidity limitation (20% – 80%)
	Unique device identifier		Do not use if package is damaged and consult instruction for use
	Responsible Person in the United Kingdom		Do not reuse
	Caution		Latex free
	Do not resterilize		Consult instruction for use
	Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number		

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BULK HYPODERMIC SYRINGE INSTRUCTION FOR USE

Product Type

Bulk Hypodermic Syringe
Luer 2 part Syringe 2 ml, 2,5 ml, 5 ml, 10 ml, 20 ml
Luer 3 part Syringe 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml
Luer Lock 3 part Syringe 1 ml, 2 ml, 2,5 ml, 3 ml, 5 ml, 10 ml, 20 ml, 50 ml, 60 ml

Syringe Type	Bulk Hypodermic Syringe
Luer – Eccentric 2 parts	
Luer – Eccentric 3 parts	
Luer – Centric 2 parts	
Luer – Centric 3 parts	
Luer Lock	

Intended Use

Bulk hypodermic syringes, made of empty, non sterile and single use plastic materials, are generally used to administer or withdraw fluids, with the intended use determined by the end user. Sterilization is performed by the end user. The product is suitable for EO sterilization. The responsibility for the sterilization process with validated parameters lies with the end manufacturer.

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Instruction for Use

The following instructions for use are intended for general use of the syringes, independent of the intended use specified by the end manufacturer.

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V*: Volume

Intended Patient Population

The patient population for the intended use of the finished product is determined by the end manufacturer.

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It is suitable for use with connection connectors according to "EN ISO 80369-7 Small – bore connectors for liquids and gases in healthcare applications – Part 7: Connectors for intravascular or hypodermic applications"

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Contraindication

It is not intended for use other than that specified by the end manufacturer.
After usage, dispose of the product as medical waste.

Cautions

- Check the product if the transport carton is damaged.
- Do not use the product after the expiration date printed on label.
- If there is a special warning for the drug administered by syringe, it should be taken into account during application.
- In applications, the capacity tolerance table of the graduated scale given in the instructions for use should be taken into account.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused, rework or re-sterilized.
- Reuse of the product may cause infections, cross-contamination and sepsis.

- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.
- The products are suitable for sterilization with EO sterilization method.

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