TR TÜBERKÜLİN ŞIRINGA KULLANIM TALİMATI

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
Tel: 0212 622 04 00

www.setmedikal.com.tr

GB TUBERCULIN SYRINGE INSTRUCTION FOR USE**Product Type**

Tuberculin Syringe with needle	Hypodermic Needle
Luer 3 pieces Syringe 1 ml	21G 26G 27G

Product Type	Tuberculin Syringe with needle – Regular Needle
Luer – Centric 3 parts	

Intended Use

Hypodermic Tuberculin Syringes are made of empty, sterile and disposable plastic materials, are used by end users for the aspiration and injection of liquids after filling. The syringes are primarily for human use. Sterile syringes are designed to be used immediately after filling and are not intended to retain medication for long period of time.

Tuberculin syringes are suitable for use with luer tip connection connectors in accordance with EN ISO 80369-7.

Sterile single use tuberculin syringes are intended to be used to administer a precise, low-dosage injection of tuberculin and/or sensitivity test allergen to a patient subcutaneously. A non-retractable needle is included, either attached or detached (usually capped for user protection). It is not prefilled with any drug.

Product Components

Name	Raw material
Barrel	PP – Polypropylene
Plunger	PP – Polypropylene
Gasket	Isoprene rubber, Latex free.
Hypodermic Needle	Stainless Steel, SS304, PP – Polypropylene
Packaging	Medical grade paper and transparent film for EtO (individually blister packed)

**Instruction for Use**

- Open the package by peeling it from the place marked by the arrow.
- Visually inspect that the content of package is complete and the package is undamaged.
- The product packaging includes the syringe and needle together.
- During priming, hold the syringe with one hand and the needle with the other.
- Make sure that the grading lines on the syringe cylinder are legible.
- Without removing the protective cap of the needle, firmly fix the syringe and needle together by pushing the needle towards the syringe tip and turning it clockwise.
- Carefully remove the needle's protective cap without damaging the needle.
- If liquid is to be injected into the body with a syringe; draw the required amount of liquid into the syringe.
- Make sure the syringe is free of air bubbles. If there is any air bubble in the syringe, be sure to remove it.
- Before injection, check if fluid is flowing in the syringe.
- Inject the liquid to be injected into the vein, intramuscularly, subcutaneously or subcutaneously, in accordance with the application procedure.
- The syringe can also be used for blood aspiration with an appropriately sized hypodermic needle.
- The electronic version of PR.A98 Instruction for use is available at www.setmedikal.com.tr

- Determine the amount of liquid to be used using the graduated scale on the syringe.
- The capacity tolerance of the graduated scale on the syringe is stated in the table below in accordance with EN ISO 7886-1 Standard.

Nominal capacity of syringe V*, (ml)	Graduated Capacity Tolerance (%)
1≤V≤3	±5

V*: Volume

Intended Patient Population

There is no limitation for the patient population.

Intended User

Healthcare personnel

Device Description Intended for Use With The Device

Suitable for use with connectors "EN ISO 80369-7 Small diameter connectors for liquids and gases in healthcare applications - Part 7: Connectors for intravenous or hypodermic applications".

Indications

Tuberculin syringes are devices used to administer medical products and liquids. The product shelf life is 5 years.

Contraindications

Tuberculin syringes with needle should not be used in more than one patient. The product can be used in a one-time sterile condition with a single patient. Cleaning or resterilization is not possible. It is not intended for use other than its stated intended use. Dispose of as medical waste after use.

Cautions

- Do not use the product if the package is damaged.
- Visually inspect that the content of package is complete and the package is undamaged.
- Do not use the product after the expiration date printed on the primary packaging.
- When the sterile barrier system (unit packaging) is opened, the needle is inserted without waiting in accordance with aseptic rules.
- After the medicine is drawn into the syringe, it should be used within a maximum of 60 minutes with the needle cap attached.
- If there is a special warning for the drug administered by syringe, it should be taken into account during application.
- In applications, the capacity tolerance table of the graduated scale given in the instructions for use should be taken into account.
- Before use, check that the hypodermic needle is not blocked and the flow is not obstructed.
- After injection or aspiration, dispose of the product as medical waste.
- The product is for single use only. It cannot be reused or re-sterilized.
- Do not touch the needle tip to prevent contamination/ needle sticks. Skin puncture with a contaminated needle may result in serious illness such as hepatitis, AIDS and known/unknown diseases.
- If the sterile packaging is not used immediately after opening, there is a risk of contamination. Product sterility cannot be maintained, do not use the product.
- Reuse of the product may cause infections, cross-contamination and sepsis.
- If any serious incident occurs, it should be reported to the manufacturer and the competent authority.


Hypodermic Needle Configurations

Needle Configurations			
Gauge	21G	26G	27G
Color	Deep Green	Brown	Medium Grey
Needle Length	5/8"	1/2"	1 1/2"
	1 1/2"		2"

Storage Conditions

- Store at ambient temperature not exceeding 45°C.
- Keep it away from direct sunlight and moist environment.
- Place and keep the box in the direction of arrow.
- Do not place more than 5 successive boxes on the top of each other.
- The shelf life of product is 5 years.

Description of Symbols Used

	Manufacturer		Fragile. Handle with care
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Move vertically
	Catalogue number		Upper limit of temperature (45°C)
	Medical device		Humidity limitation (20% – 80%)
	Unique device identifier		Do not use if package is damaged and consult instruction for use
	Responsible Person in the United Kingdom		Do not reuse
	Sterilized using ethylene oxide		Latex free
	Do not resterilize		Non-pyrogenic
	Single sterile barrier system		Professional Use Only
	Single sterile barrier system with protective packaging outside		Consult instruction for use or consult electronic instruction for use

	Caution		Marking for meeting the requirements of UK Medical Device Regulation The number displayed on the bottom is the Approved Body Number
	Marking for meeting the requirements of all relevant European Medical Device Directives/European Legislation The number displayed on the bottom is the Notified Body Number		

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