





set[®] diabetes

Daily diabetes management relies on tools that feel trustworthy and familiar.

Our diabetes care range focuses on comfort, precision, and ease of use, offering insulin syringes and needles

that support users in their everyday routines.

These products are developed to deliver consistent dosing while fitting seamlessly into the rhythm of self-care.




We care.

Diabetes Instruction

Our DIABETES offerings

For syringes, proper dosage accuracy and a needle tip that meets standard requirements are of great importance.



Various medical devices are used during the diabetes management process. These include syringes designed for accurate dosing and insulin pens used for insulin administration.



Thin walled body and specific angled sharp tip that minimizes pain



For insulin pens, a pen needle is attached to the tip of the pen. This product is essential for subcutaneous injection of insulin and helps make the application process easier.



set[®] diabetes

Insulin Syringes



We care.

Set® Insulin

Features

- Luer Slip tip, compatible with U-100 insulin
- Precise and smooth plunger movement with medical-grade silicone lubrication
- Excellent sealing performance with latex-free isoprene rubber gasket
- No air leakage through gasket or needle junction
- Clear and legible graduation (0.01 ml / 0.02 ml depending on model)
- Fixed needle (26G / 29G / 30G / 31G) or 6% Luer tip depending on syringe type
- Transparent barrel for easy insulin level visibility
- Designed for safe, single-handed subcutaneous insulin administration
- Durable structure compliant with intended injection pressure


Packaging

- Sterile or non-sterile
- With single or double cap
- Without or with needle. (bi-packed)
- Without or with aspiration or filter needle

Conformance & Compliance

- ISO 8537 – ISO 80369-7 – ISO 9626
- EU 2017/745 (MDR)
- FDA 510(k) Cleared
- EO sterilized
- Shelf life: 5 years

Available Sizes & Types



26G x 1/2"
0.45 x 13 mm




Double Cap
Insulin Syringe



29G x 1/2"
0.33 x 13 mm




Single Cap
Insulin Syringe



30G x 5/16"
0.30 x 8 mm



31G x 5/16"
0,25x12mm



Product Table

Product Description	Syringe Volume	Units / Box	Units / Carton
0,3ml 3-Part 29Gx1/2	0.3 ml	100/450 pcs	2400/2700 pcs
0,3ml 3-Part 30Gx5/16	0.3 ml	100/450 pcs	2400/2700 pcs
0,3ml 3-Part 31Gx5/16	0.3 ml	100/450 pcs	2400/2700 pcs
0,5ml 3-Part 29Gx1/2	0.5 ml	100/450 pcs	2400/2700 pcs
0,5ml 3-Part 30Gx5/16	0.5 ml	100/450 pcs	2400/2700 pcs
0,5ml 3-Part 31Gx5/16	0.5 ml	100/450 pcs	2400/2700 pcs
01ml 3-Part 29Gx1/2	1 ml	100/450 pcs	2400/2700 pcs
01ml 3-Part 30Gx5/16	1 ml	100/450 pcs	2400/2700 pcs
01ml 3-Part 30Gx1/2	1 ml	100/450 pcs	2400/2700 pcs
01ml 3-Part 26Gx1/2	1 ml	100/250 pcs	2400/3000 pcs

set[®] diabetes

Pen Needle



We care.



Set® Lance


Features

- Manufactured with SS304 stainless-steel needle
- Atraumatic, triple-bevel tip design for smooth and comfortable penetration
- Silicone lubrication ensuring consistent, low-friction needle entry
- Polypropylene hub, protective cap and needle shield
- Suitable for use with all standard insulin pens
- Available in multiple gauges and lengths for patient-specific dosing (30G–32G, 4–8 mm)
- Designed for safe, single-use insulin injection

Conformance & Compliance

- 93/42/EEC (MDD)
- ISO 7864 - ISO 11608-2
- EO sterilized
- Shelf life: 5 years

Available Sizes & Types



Thin walled body and specific angled sharp tip that minimizes pain



Product Table

Size	Needle Gauge (G)	Units / Box	Units / Carton
30G x 8mm	30	100 pcs	6000 pcs
31G x 4mm	31	100 pcs	6000 pcs
31G x 5mm	31	100 pcs	6000 pcs
31G x 6mm	31	100 pcs	6000 pcs
31G x 8mm	31	100 pcs	6000 pcs
32G x 4mm	32	100 pcs	6000 pcs
32G x 6mm	32	100 pcs	6000 pcs
32G x 8mm	32	100 pcs	6000 pcs