



set[®] catheters

Catheterization is central to many diagnostic and therapeutic procedures, requiring devices that work smoothly and predictably.

Our catheter portfolio includes a variety of solutions for access, drainage, and fluid management, all designed to support efficient workflows and dependable performance in daily practice.



We care.

Catheters Instruction

Our CATHETER offerings

Procedures can be performed using different vacuum catheters for different purposes. Some catheters are vacuum-controlled.

They can be available in different types and forms, ranging from airway clearance to rectal catheters or urinary bladder catheters.



Rectal Catheters - Tiemann Catheter - Nelaton Catheters - Suction Catheters





Set® Cath Respiratory Care

Suction Catheters



Features

- Smooth surface enabling easy and controlled insertion
- Optimal stiffness preventing kinking during use
- Atraumatic, rounded open distal tip minimizing tissue trauma
- Two side holes designed to prevent tissue aspiration (cross-oriented holes available in specific models)
- Available tip designs: Bevel Tip and Straight / Ideal Tip
- Color-coded funnel connector for quick size identification
- Optional vacuum control connector providing effective suction regulation
- Universal connector for all models
- Standard catheter length: 50 cm

Conformance & Compliance

- 93/42/EEC (MDD)
- ISO 8836

Available Sizes & Types



Product Table

Size	Ideal Tip	Beveled Tip	Units Per Box (Vacuum Control)	Units Per Box (Regular)
● ch/fr 10	●	●	100pcs	150pcs
● ch/fr 12	●	●	100pcs	150pcs
● ch/fr 14	●	●	100pcs	150pcs
● ch/fr 16	●	●	100pcs	150pcs
● ch/fr 18	●	●	100pcs	150pcs



Set® Cath Urology Care

Tiemann Catheters

Features


- Designed for effective catheterization in obstructed or narrow passages
- Smooth surface enabling easy and controlled insertion
- Optimal stiffness preventing kinking during advancement
- Atraumatic, rounded angled closed distal tip minimizing tissue irritation
- Cross-positioned double side holes helping maintain uninterrupted flow
- Color-coded connector for quick size identification
- Standard catheter length: 40

Conformance & Compliance


- 93/42/EEC (MDD)




Available Sizes & Types




ch/fr 10




ch/fr 12



ch/fr 14



ch/fr 16



ch/fr 18

Product Table

Size	Units Per Box	Units Per Case
● ch/fr 10	100pcs	1000pcs
● ch/fr 12	100pcs	1000pcs
● ch/fr 14	100pcs	1000pcs
● ch/fr 16	100pcs	1000pcs
● ch/fr 18	100pcs	1000pcs



Set® Cath Gastrointestinal Care


Rectal Catheters

Features


- Smooth surface enabling easy and controlled insertion
- Optimal stiffness preventing kinking during use
- Atraumatic, rounded closed distal tip minimizing tissue irritation
- Cross-positioned double side holes designed to improve fluid drainage
- Color-coded connector for quick size identification
- Standard catheter length: 40 cm

Conformance & Compliance


- 93/42/EEC (MDD)




Available Sizes & Types




ch/fr 18




ch/fr 20



ch/fr 22



ch/fr 24



ch/fr 26

Product Table

Size	Units Per Box	Units Per Case
ch/fr 18	75pcs	750pcs
ch/fr 20	75pcs	750pcs
ch/fr 22	75pcs	750pcs
ch/fr 24	75pcs	750pcs
ch/fr 26	75pcs	750pcs

